

Recombinant Human Cdc26 / Apc12 protein ab108123

1 图像

描述	
产品名称	重组人Cdc26 / Apc12蛋白
纯度	> 90 % SDS-PAGE. ab108123 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q8NHZ8</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMLRRKPTRLCLKDDIEEFE NIRKDLETRK KQKEDVEVVGGSDGEGAIGLSSDPKSREQMINDRIGYKPQPK PNNRSSQF GSLEF
预测分子量	12 kDa including tags
氨基酸	1 to 85
标签	His tag N-Terminus
技术指标	
Our Abpromise guarantee covers the use of ab108123 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


pH: 8.00

Constituents: 0.316% Tris HCl, 40% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

相关性	The protein encoded by this gene is highly similar to <i>Saccharomyces cerevisiae</i> Cdc26, a component of cell cycle anaphase-promoting complex (APC). APC is composed of a group of highly conserved proteins and functions as a cell cycle-regulated ubiquitin-protein ligase. APC thus is responsible for the cell cycle regulated proteolysis of various proteins.
细胞定位	Nuclear

图片



15% SDS-PAGE analysis of 3µg ab108123.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors