

Recombinant Human CEACAM1 protein ab151901

描述		
产品名称	重组人CEACAM1蛋白	
纯度	> 95 % SDS-PAGE. Purity is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.	
内毒素水平	< 1.000 Eu/μg	
表达系统	HEK 293 cells	
Accession	<u>P13688</u>	
蛋白长度	Protein fragment	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	QLTTESMPFNVAEGKEVLLLLVHNLPQQLFGYSWYKGERVDGN RQIVGYAI GTQQATPGPANSGRETIYPNASLLIQNVTQNDTGFYTLQVIK SDLVNEEA TGQFHVYPELPKPSISSNNSNPVEDKDAVAFTCEPETQDTTY LWWINNQS LPVSPRLQLSNGNRTLTLSSVTRNDTGPYECEIQNPVSANRS DPVTLNVT YGPDTPTISPSDTYRPGANLSLSCYAASNPPAQYSWLINGT FQQSTQEL FIPNITVNNSGSYCHANNSVTGCNRTTVKTIIVTELSPVVA KPQIKASK TTVTGDKDSVNLTCSTNDTGISIRWFFKNQSLPSSERMKLSQ GNTTLSIN PVKREDAGTYWCEVFNPISKQSDPIMLNVNYNALPQENGLS PGVDHHHH HH	
预测分子量	44 kDa including tags	
氨基酸	35 to 428	
标签	His tag C-Terminus	

技术指标

Our **Abpromise guarantee** covers the use of **ab151901** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
	HPLC
形式	Lyophilized

制备和贮存

稳定性和存储	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.
	pH: 7.20 Constituents: 99% Phosphate Buffer, 0.88% Sodium chloride
复溶	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in 1X PBS. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. For long term storage aliquot and store at < -20°C.

常规信息

序列相似性	Belongs to the immunoglobulin superfamily. CEA family.
	Contains 3 Ig-like C2-type (immunoglobulin-like) domains.
	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Secreted; Cell membrane and Cell membrane. Localizes to sites of cell-cell contact.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors