

Recombinant Human Clathrin light chain B protein ab130049

1 图像

描述	
产品名称	重组人Clathrin light chain B蛋白
纯度	> 85 % SDS-PAGE. ab130049 was purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>P09497</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MGSMADDFGF FSSSESGAPE AAEDPAAAF LAQGESEIAG IENDEGFGAP AGSHAAPAQP GPTSGAGSED MGTTVNGDVF QEANGPADGY AAIAQADRLT QEPESIRKWR EEQRKRLQEL DAASKVTEQE WREKAKKDLE EWNQRQSEQV EKNKINNRAS EEAFVKESKE ETPGTEWEKV AQLCDFNPKS SKQCKDVSRL RSVLMSLKQT PLSR
预测分子量	26 kDa
氨基酸	1 to 211
标签	His tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab130049 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Mass Spectrometry SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.004% PMSF, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride


常规信息

功能 Clathrin is the major protein of the polyhedral coat of coated pits and vesicles.

序列相似性 Belongs to the clathrin light chain family.

细胞定位 Cytoplasmic vesicle membrane. Membrane > coated pit. Cytoplasmic face of coated pits and vesicles.

图片



15% SDS-PAGE analysis of 3 µg of Human Clathrin light chain B protein (ab130049)

Note: Molecular weight on SDS-PAGE will appear higher than predicted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors