

Recombinant Human DNTTIP1/TDIF1 protein ab171701

1 图像

描述

产品名称	重组人DNTTIP1/TDIF1蛋白
纯度	> 85 % SDS-PAGE. ab171701 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q9H147</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<pre> MGSSHHHHHH SSGLVPRGSH MGSMGATGDA EQPRGPSGAE RGGLELG DAG AAGQLVLTNP WNIMIKHRQV QRRGRRSQMT TSFTDPAISM DLLRAVLQPS INEEIQT VFN KYMKFFQKAA LNVRDNV GEE VDAEQLIQEA CRSCLEQAKL LFSDG EKVIP RLTHELPGIK RGRQAEEECA HRGSPLPKKR KGRPPGHILS SDRAAAGMVW KPKSCEPIRR EGPKWDPARL NESTTFVLGS RANKALGMGG TRGRIYIKHP HLFKYAADPQ DKHWLAEQHH MRATGGKMAY LLIEEDIRDL AASDDYRGCL DLKLEELKSF VLPSWMVEKM RKYMETLRTE NEHRAVEAPP QT </pre>
预测分子量	39 kDa including tags
氨基酸	1 to 329
标签	His tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab171701** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
	Mass Spectrometry
质谱法	MALDI-TOF

形式	Liquid
补充说明	This product was previously labelled as DNTTIP1


制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.02% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.88% Sodium chloride
---------------	---

常规信息

功能	Shown to enhance TdT activity, in vitro.
细胞定位	Nucleus.

图片



15% SDS-PAGE analysis of ab171701 (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors