

Recombinant Human Galectin 13/PLAC8 protein ab124309

1 图像

描述

产品名称	重组人Galectin 13/PLAC8蛋白
纯度	> 90 % SDS-PAGE. ab124309 was purified using Ni-NTA chromatography. The product is filtered (0.4 µm) but not sterile, hence needs to be filtered by an appropriate sterile filter before using in cell culture.
表达系统	Escherichia coli
Accession	<u>Q9UHV8</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<pre> MKHHHHHHAS MSSLPVPYKL PVSLSVGSCV IIKGTPIHFS INDPQLQVDF YTDMEDSDI AFRFRVHFGN HVVMNRREFG IWMLEETDY VPFEDGKQFE LCIYVHYNEY EIKVNGIRIY GFVHRIPPSF VKMVQVSRDI SLTSVCVCN </pre>
预测分子量	17 kDa including tags
氨基酸	1 to 139
标签	His tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab124309** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot ELISA SDS-PAGE
形式	Lyophilized
补充说明	This product was previously labelled as Galectin 13

制备和贮存

稳定性和存储

Shipped at 4°C. Store at -80°C.

Constituents: 99% Phosphate Buffer, 0.44% Sodium chloride

复溶

Add deionized water to prepare a working stock solution of 0.5 mg/mL and let the lyophilized pellet dissolve completely. Avoid storage of lyophilized product at room temperature. Aliquot reconstituted protein to avoid repeated freezing/thawing cycles and store at -80°C for long term. Reconstituted protein can be stored at 4°C for a week.

常规信息

功能

Has lysophospholipase activity.


组织特异性

Detected in adult and fetal spleen, fetal kidney, adult urinary bladder and placenta.

序列相似性

Contains 1 galectin domain.

图片



14% SDS-PAGE separation of Human Galectin 13/PLAC8 protein (ab124309)

Lane 1: reduced and boiled sample (5 µg)

Lane 2: non-reduced and non-boiled sample (5 µg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors