

Recombinant Human GATA2 protein ab134866

1 References

描述

| | |
|------------------|--|
| 产品名称 | 重组人GATA2蛋白 |
| 纯度 | > 90 % SDS-PAGE. |
| 表达系统 | Escherichia coli |
| Accession | <u>P23769</u> |
| 蛋白长度 | Full length protein |
| 无动物成分 | No |
| 性质 | Recombinant |
| 种属 | Human |
| 序列 | (29 amino acid tag)EVAPEQPRWMAHPAVLNAQHPDSSHHPGLAHNYMEPAQL LPPDEV DV FFNHLD SQGNPYANPAHARARVSYSPA HARLTGGQMCRPHL LHSPGLPW LDGGKAALSAAAAHHHNPWTVSPFSKTPLHPSAAGGGGGLS VYPGAGGG SGGGSGSSVASLTPTAAHSGSHLFGFPPTPPKEVSPDPSTTG AASPASSS AGGSAARGEDKDGVKYQVSLTESMKMESGSPLRPGLATMG TQ PATHHPIP TYP SYV PAAAH DYSSGLFHPGGFLGGPASSFTPKQRSKARSC SEGRECVN CGATATPLWRRDGTGHYLCNACGLYHKMNGQNRPLIKPKRRL SAARRAGT CCANCQTTTTTLWRRNANGDPVCNACGLYYKLHNVRPLTMK KEGIQTRN RKMSNKSKKSKGAECFEELSKCMQEKSSPFSAALAGHMAP VGHLP PFS HSGHILPTPTPIHPSSSL SFGHPHPSSMVTAMGLEESGGGGS PGRRRRRR RRRRR |
| 预测分子量 | 50 kDa |
| 氨基酸 | 2 to 480 |

技术指标

Our **Abpromise guarantee** covers the use of **ab134866** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----|----------|
| 应用 | SDS-PAGE |
| 形式 | Liquid |

制备和贮存

| | |
|--------|---|
| 稳定性和存储 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. pH: 8.00 Constituents: Potassium chloride, 0.05% DTT, 0.32% Tris HCl, 0.02% EDTA, Glycerol, Sodium chloride, 3.4% DL-Arginine |
|--------|---|

常规信息

| | |
|-------|--|
| 功能 | Transcriptional activator which regulates endothelin-1 gene expression in endothelial cells. Binds to the consensus sequence 5'-AGATAG-3'. |
| 组织特异性 | Endothelial cells. |
| 序列相似性 | Contains 2 GATA-type zinc fingers. |
| 细胞定位 | Nucleus. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors