

Recombinant Human GIP protein ab151376

1 References

描述	
产品名称	重组人GIP蛋白
纯度	> 95 % SDS-PAGE.
内毒素水平	< 1.000 Eu/µg
表达系统	Cell Culture
Accession	<u>Q9BXJ1</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	RVPHVQGEQQEWEGTEELPSPPDHAERAEQHEKYRPSQDQG LPASRCLR CCDPGTSMYP ATAVPQINITILKGEKGDRGDRGLQGKYGKTGSAGARGHTGP KGQKGSMD APGERCKSHY AAFSVGRKKPMHSNHYYQTVIFDTEFVNLYDHFNMFTGKFYC YVPGLYFF SLNVHTWNQK ETYLHIMKNEEEVVILFAQVGDRSIMQSQSLMLELREQDQVW VRLYKGER ENAIFSEELD TYITFSGYLVKHATEP
预测分子量	29 kDa
氨基酸	26 to 281
标签	His tag C-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab151376 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE HPLC
形式	Lyophilized
补充说明	Previously labelled as CTRP1.

制备和贮存

稳定性和存储

Shipped at 4°C. Store at -20°C.

pH: 7.2

复溶

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles

常规信息

序列相似性

Contains 1 C1q domain.

Contains 1 collagen-like domain.

细胞定位

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors