

Recombinant human ICK protein ab101774

4 图像

描述

产品名称	重组人ICK蛋白
生物活性	The Specific activity of ab101774 was determined to be 4.8 nmol/min/mg.
纯度	> 75 % SDS-PAGE. Purity was determined to be >75% by densitometry. Affinity purified.
表达系统	Baculovirus infected insect cells
Accession	Q9UPZ9
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	98 kDa including tags
氨基酸	1 to 632

技术指标

Our [Abpromise guarantee](#) covers the use of **ab101774** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot Functional Studies SDS-PAGE
形式	Liquid
补充说明	ab64311 (Myelin Basic Protein protein) can be utilized as a substrate for assessing kinase activity

制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.307% Glutathione, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
--------	--


This product is an active protein and may elicit a biological response in vivo, handle with caution.

常规信息


相关性 ICK may play a key role in the development of multiple organ systems and particularly in cardiac development.

细胞定位 Cytoplasmic and Nuclear


图片



The specific activity of ICK (ab101774) was determined to be 4.8 nmol/min/mg as per activity assay protocol




SDS PAGE analysis of ab101774



SDS-PAGE analysis of ab101774.

SDS-PAGE - Recombinant human ICK protein
(ab101774)



The Specific activity of ab101774 was determined to be 4.8 nmol/min/mg.

Functional Studies - Recombinant human ICK
protein (ab101774)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors