

Recombinant Human IGSF10 protein ab166199

1 图像


描述	
产品名称	重组人IGSF10蛋白
表达系统	Wheat germ
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	SAFISPQGFMAPFGSLTLNMTDQSGNEANMVCSIQKPSRTSP IAFTEEND YIVLNTSFSTFLVCNIDYGHIQPVWQILALYSDSPLILERSH LLSETPQ
氨基酸	301 to 399
标签	GST tag N-Terminus

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab166199</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	ELISA Western blot
形式	Liquid
补充说明	

制备和贮存	
稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl

功能	May be involved in the maintenance of osteochondroprogenitor cells pool.
序列相似性	Contains 12 Ig-like C2-type (immunoglobulin-like) domains. Contains 6 LRR (leucine-rich) repeats. Contains 1 LRRCT domain. Contains 1 LRRNT domain.
细胞定位	Secreted.

## 图片



ab166199 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors