

Recombinant Human IL-18BP ab193078

描述

产品名称	重组人IL-18BP
纯度	> 95 % SDS-PAGE. Purity greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
内毒素水平	< 1.000 Eu/μg
表达系统	HEK 293 cells
Accession	<u>O95998</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	TPVSQTTTAATASVRSTKDPSPQPPVFPAKQCPALEVTWP EVEVPLNG TSLSLSCVACSRFPNFSILYWLGNFSFIEHLPGRLEWEGSTSR RGSTGTQL CKALVLEQLTPALHSTNFCVLDPEQVVQRHVLAQLWAGL RATLPPTQ EALPSSHSSPQQGVDDIEGRMDEPKSCDKTHTCPPCPAPEL LGGPSVFL FPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEV HNAKTKPR EEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEK TISKAKGQ PREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESN GQPENNYK TTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHN HYTQKSLS LSPGK
预测分子量	45 kDa including tags
氨基酸	31 to 194
标签	Fc tag C-Terminus
额外的序列信息	This product is the mature full length protein from aa 31 to 194. The signal peptide is not included. Fused with an Fc tag at the C terminus.
描述	重组人IL-18BP蛋白

技术指标

Our **Abpromise guarantee** covers the use of **ab193078** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
	HPLC
形式	Lyophilized

制备和贮存

稳定性和存储 Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituent: 100% PBS

Lyophilized from a 0.2 µM filtered solution.

复溶 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

常规信息

功能	Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.
组织特异性	Strongly expressed in heart, lung, placenta and spleen.
序列相似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors