

Recombinant Human Inhibin beta A protein ab53506

1 图像

描述

产品名称	重组人Inhibin beta A蛋白
纯度	> 90 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<u>P08476</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
氨基酸	311 to 426
标签	His tag N-Terminus
额外的序列信息	The sequence is the mature full length without the signal peptide and propeptide. Hexahistidine tag at N-terminus.

技术指标

Our **Abpromise guarantee** covers the use of **ab53506** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
--------	---

常规信息

功能	Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival,
----	---

embryonic axial development or bone growth, depending on their subunit composition. Inhibins appear to oppose the functions of activins.


序列相似性

Belongs to the TGF-beta family.

细胞定位

Secreted.

图片



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors