

Recombinant Human Interferon alpha 8 protein (His tag) ab276257

1 图像

描述

产品名称	重组人Interferon alpha 8蛋白(His tag)		
纯度	> 90 % SDS-PAGE.		
内毒素水平	< 1.000 Eu/μg		
表达系统	HEK 293 cells		
Accession	<u>P32881</u>		
蛋白长度	Full length protein		
无动物成分	No		
性质	Recombinant		
种属	Human		
序列	MALTFYLLVA	LVVLSYKSFS	SLGCDLPQTH
	SLGNRRALIL	LAQMRRISPF	SCLKDRHDFE
	FPQEEFDDKQ	FQKAQAISVL	HEMIQQTfNL
	FSTKDSSAAL	DETLLEDFYI	ELDQQLNDLE
	SCVMQEVGVI	ESPLMYEDSI	LAVRKYFQRI
	TLYLTEKKYS	SCAWEVVRAE	IMRSFSL SIN LQKRLKSKE
预测分子量	21 kDa including tags		
分子量信息	The recombinant human IFNA8 comprises 177 amino acids and has a predicted molecular mass of 20.9 kDa. The apparent molecular mass of the protein is approximately 28 kDa in SDS-PAGE under reducing conditions.		
氨基酸	1 to 189		
标签	His tag C-Terminus		
额外的序列信息	Predicted N-terminus: Cys 24.		

技术指标

Our **Abpromise guarantee** covers the use of **ab276257** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 SDS-PAGE

形式 Lyophilized

制备和贮存

稳定性和存储 Shipped at Room Temperature. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituent: 100% PBS


复溶 This information is lot specific. Please contact our technical Support team for details.

常规信息

相关性 Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

细胞定位 Secreted

图片



SDS-PAGE analysis of ab276257

SDS-PAGE - Recombinant Human Interferon alpha 8 protein (ab276257)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors