

Recombinant Human Isoleucyl tRNA synthetase protein

ab152465

1 图像

描述	
产品名称	重组人Isoleucyl tRNA synthetase蛋白
表达系统	Wheat germ
Accession	<u>P41252</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	QYINLQLLNAKPQECLMGTVGTLLENPLGQNGLTHQGLLYE AAKVFGLR SRKLLKFLNETQTQEITEDIPVKTLNMKTVYVSVLPTTADF
预测分子量	36 kDa including tags
氨基酸	1172 to 1262


技术指标	
Our <u>Abpromise guarantee</u> covers the use of ab152465 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot SDS-PAGE ELISA
形式	Liquid
补充说明	

制备和贮存	
稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl

常规信息

相关性	Isoleucyl tRNA synthetase is a component of multi ARS complex and catalyzes the ligation of isoleucine with isoleucyl tRNA.
细胞定位	Cytoplasmic

图片



12.5% SDS-PAGE analysis of ab152465 stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors