

Recombinant Human MESH1 protein ab113583

1 图像

| | |
|-----------|---|
| 描述 | |
| 产品名称 | 重组人MESH1蛋白 |
| 纯度 | > 90 % SDS-PAGE. ab113583 is purified using conventional chromatography techniques: anion exchange chromatography followed by gel filtration chromatography with 20mM Tris pH 7.5, 2mM EDTA. No detergents such as urea, Triton or Tween were used to purify the proteins. They purified the proteins by using |
| 表达系统 | Escherichia coli |
| Accession | <u>Q8N4P3</u> |
| 蛋白长度 | Full length protein |
| 无动物成分 | No |
| 性质 | Recombinant |
| 种属 | Human |
| 序列 | MGSSHHHHHHSSGLVPRGSHLDEVELHFGAQVRRLLVEEVT DDKTLPKLER KRLQVEQAPHSSPGAKLVKLADKLYNLRDLNRCTPEVKIQMG SEAAQLLE AADFAARKHRQRRKDPEGTPYINHPIGVARILTHEAGITDI VVLQAALL HDTVEDTDTT |
| 预测分子量 | 18 kDa including tags |
| 氨基酸 | 1 to 140 |
| 标签 | His tag N-Terminus |

| | |
|--|-------------------------------|
| 技术指标 | |
| Our Abpromise guarantee covers the use of ab113583 in the following tested applications. | |
| The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user. | |
| 应用 | SDS-PAGE Mass Spectrometry |
| 质谱法 | MALDI-TOF |
| 形式 | Liquid |
| 补充说明 | Previously labelled as HDDC3. |

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 40% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

常规信息

功能


ppGpp hydrolyzing enzyme involved in starvation response.

序列相似性

Belongs to the MESH1 family.

Contains 1 HD domain.

图片



15% SDS-PAGE showing ab113583 at approximately 17.9 kDa (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors