

Recombinant Human Miz1 protein ab151354

描述		
产品名称	重组人Miz1蛋白	
纯度	> 95 % SDS-PAGE.	
内毒素水平	< 1.000 Eu/μg	
表达系统	Escherichia coli	
Accession	<u>Q13105</u>	
蛋白长度	Protein fragment	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	MGSSHHHHHHSSGLVPRGSHMDFPQHSQHVLEQLNQQRQLGL LCDCTFVV DGVHFKAHKAVLAACSEYFKMLFVDQKDVVHLDISNAAGLGQ VLEFMYTA KLSLSPENVDDVLAVATFLQMQDIITACHALKSLAEPATSPG GNAEALAT EGGDKRAKEEKVATSTLSRLEQAGRSTPIGPSRDLKEERGGQ AQSAASGA EQTEKADA	
预测分子量	22 kDa including tags	
氨基酸	1 to 188	
标签	His tag N-Terminus	

技术指标		
Our Abpromise guarantee covers the use of ab151354 in the following tested applications.		
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.		
应用	SDS-PAGE	
形式	Lyophilized	

制备和贮存		
稳定性和存储	Shipped at 4°C. Store at -80°C.	

pH: 7.4

Constituents: 0.88% Sodium chloride, 99% Phosphate Buffer

复溶

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in 1X PBS. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

常规信息

功能

May function as a housekeeping DNA-binding protein that regulates the expression of specific genes. Has been shown to bind to the promoters of adenovirus major late protein and cyclin D1 and activate transcription. Also has potent growth arrest activity, probably through inhibition of cell cycle progression. Required for early embryonic development during gastrulation.

序列相似性

Belongs to the krueppel C2H2-type zinc-finger protein family.
Contains 1 BTB (POZ) domain.
Contains 13 C2H2-type zinc fingers.

细胞定位

Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors