

Recombinant Human OXSR1 protein ab125579

1 图像

描述

产品名称	重组人OXSR1蛋白
纯度	> 95 % Densitometry. Purity was determined to be >95% by densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
Accession	<u>O95747</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	86 kDa including tags
氨基酸	1 to 527

技术指标

Our **Abpromise guarantee** covers the use of **ab125579** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot SDS-PAGE
形式	Liquid

制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.29% Sodium chloride
--------	---

常规信息

功能	Regulates downstream kinases in response to environmental stress. May also have a function in
----	---

regulating the actin cytoskeleton.

组织特异性

Ubiquitously expressed in all tissue examined.


序列相似性

Belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. STE20 subfamily. Contains 1 protein kinase domain.

翻译后修饰

Phosphorylated upon DNA damage, probably by ATM or ATR.

图片



SDS-PAGE analysis of ab125579.

SDS-PAGE - Recombinant Human OXSR1 protein (ab125579)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors