

Recombinant Human PSME2 protein ab117228

1 图像

描述		
产品名称	重组人PSME2蛋白	
纯度	> 95 % SDS-PAGE. ab117228 is purified using conventional chromatography.	
表达系统	Escherichia coli	
Accession	<u>Q9UL46</u>	
蛋白长度	Full length protein	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	MGSSHHHHHHSSGLVPRGSHMAKPCGVRLSGEARKQVEVFRQ NLFQEAE FLYRFLPQKIIYLNQLLQEDSLNVADLTSLRAPLDIIPDPP PKDDEMET DKQEKKEVPKCGFLPGNEKVL SLLALVKPEVWTLKEKCILVI TWIQLHIP KIEDGNDFGVAIQEKVLERNVAVKTKVEAFQTTISKYFSERG DAVAKASK ETHVMDYRALVHERDEAAYGELRAMVLDLRAFYAELYHIIS NLEKIVNP KGEKPSMY	
预测分子量	30 kDa including tags	
氨基酸	1 to 239	
标签	His tag N-Terminus	

技术指标	
Our Abpromise guarantee covers the use of ab117228 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Mass Spectrometry SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.002% PMSF, 0.02% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息


功能

Implicated in immunoproteasome assembly and required for efficient antigen processing. The PA28 activator complex enhances the generation of class I binding peptides by altering the cleavage pattern of the proteasome.

序列相似性

Belongs to the PA28 family.

图片



SDS-PAGE analysis of ab117228

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors