

Recombinant Human RELM beta protein (Animal Free) ab218075

描述

产品名称	重组人RELM beta蛋白(Animal Free)
纯度	> 95 % SDS-PAGE.
内毒素水平	< 0.500 Eu/μg
表达系统	Escherichia coli
Accession	<u>Q9BQ08</u>
蛋白长度	Full length protein
无动物成分	Yes
性质	Recombinant
种属	Human
序列	MQCSLDSVMD KIKDVLNSL EYSPSPISKK LSCASVKSQG RPSSCPAGMA VTGCACGYGC GSWDVQLETT CHCQCSVVDW TTARCCHLT
预测分子量	19 kDa
氨基酸	24 to 111
额外的序列信息	This product is for the mature full length protein. The signal peptide is not included.

技术指标

Our **Abpromise guarantee** covers the use of **ab218075** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Lyophilized

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. For long term storage it is recommended to add a carrier protein on reconstitution (0.1% HSA or BSA). Working aliquots stored with a carrier protein are stable for at least 3 months at -20°C to -80°C.. Constituent: 0.1% Trifluoroacetic acid
--------	--

Lyophilized from a sterile filtered aqueous solution.

复溶

Sterile 0.1% Trifluoroacetic Acid (TFA) at 0.1 mg/mL. Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution.

常规信息

功能	Probable hormone.
组织特异性	Expressed only in the gastrointestinal tract, particularly the colon.
序列相似性	Belongs to the resistin/FIZZ family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors