

# Recombinant Human RNF114 protein (denatured) ab177670

## 1 图像

### 描述

产品名称	重组人RNF114蛋白(denatured)
纯度	> 90 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<u><a href="#">Q9Y508</a></u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SGLVPRGSH MGSMAAQQRD CGGAAQLAGP AAEADPLGRF TCPVCLEVYE KPVQVPCGHV FCSACLQECL KPKKPVCGVC RSALAPGVRA VELERQIEST ETSCHGCRKN FFLSKIRSHV ATCSKYQNYI MEGVKATIKD ASLQPRNVPN RYTFPCPYCP EKNFDQGLV EHCKLFHSTD TKSVVCPICA SMPWGDPNYR SANFREHIQR RHRFSYDTFV DYDVDEEDMM NQVLQRSIID Q
预测分子量	28 kDa including tags
氨基酸	1 to 228
标签	His tag N-Terminus
额外的序列信息	(NP_061153)

### 技术指标

Our **Abpromise guarantee** covers the use of **ab177670** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid
补充说明	Previously labelled as ZNF313.

## 制备和贮存

---

### 稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 2.4% Urea, 10% Glycerol (glycerin, glycerine)

## 常规信息

---

### 功能

May play a role in spermatogenesis.

### 组织特异性

Expressed in numerous tissues, including skin, CD4 lymphocytes and dendritic cells. Highest levels in testis.


### 序列相似性

Contains 1 RING-type zinc finger.

---

## 图片

---



15% SDS-PAGE analysis of ab177670 (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors