

Recombinant Human RPS4X/SCAR protein ab174424

1 图像

描述

产品名称	重组人RPS4X/SCAR蛋白
纯度	> 95 % SDS-PAGE. ab174424 is purified by using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA.
表达系统	Escherichia coli
Accession	<u>NP_000998</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<pre> MGSSHHHHHH SSGLVPRGSH MGSMARGPKK HLKRVAAPKH WMLDKLTGVF APRPSTGPHK LRECLPLIIF LRNRLKYALT GDEVKKICMQ RFIKIDGKVR TDITYPAGFM DVISIDKTGE NFRLIYDTKG RFAVHRITPE EAKYKLCKVR KIFVGTGKIP HLVTHDARTI RYPDPLIKVN DTIQIDLETG KITDFIKFDT GNLCMVTGGA NLGRIGVITN RERHPGSFDV VHVKDANGNS FATRLSNIFV IGKGNKPWIS LPRGKGIRLT IAEERDKRLA AKQSSG </pre>
预测分子量	32 kDa including tags
氨基酸	1 to 263
标签	His tag N-Terminus
额外的序列信息	X isoform.

技术指标

Our **Abpromise guarantee** covers the use of **ab174424** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Mass Spectrometry
	SDS-PAGE

质谱法	MALDI-TOF
形式	Liquid
补充说明	This product was previously labelled as RPS4X

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 0.88% Sodium chloride, 40% Glycerol (glycerin, glycerine), 0.02% DTT

常规信息

序列相似性


Belongs to the ribosomal protein S4e family.

Contains 1 S4 RNA-binding domain.

细胞定位

Cytoplasm. Localized in cytoplasmic mRNP granules containing untranslated mRNAs.

图片



SDS-PAGE analysis of 3 µg ab174424 on 15% SDS-PAGE gel.

SDS-PAGE - Recombinant Human RPS4X/SCAR protein (ab174424)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors