

Recombinant Human Serpin A5 protein ab89362

1 图像

描述

| | |
|-------|---------------------------------------------------------------------------------------|
| 产品名称 | 重组人Serpin A5蛋白 |
| 纯度 | > 90 % SDS-PAGE. ab89362 is purified using conventional chromatography techniques. |
| 表达系统 | Escherichia coli |
| 蛋白长度 | Full length protein |
| 无动物成分 | No |
| 性质 | Recombinant |
| 种属 | Human |

技术指标

Our **Abpromise guarantee** covers the use of **ab89362** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|------|---------------------------------------------|
| 应用 | SDS-PAGE |
| 形式 | Liquid |
| 补充说明 | Previously labelled as Protein C inhibitor. |


制备和贮存

| | |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 稳定性和存储 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Store in the dark. Avoid repeated freeze / thaw cycles. pH: 6 Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine) |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

常规信息

| | |
|-------|-----------------------------------------------------------------|
| 功能 | Inhibits activated protein C as well as plasminogen activators. |
| 组织特异性 | Expressed by the liver and secreted in plasma. |
| 序列相似性 | Belongs to the serpin family. |
| 细胞定位 | Secreted. |

图片



SDS-PAGE - Recombinant Human Serpin A5
protein (ab89362)

3ug by SDS-PAGE under reducing condition and visualized by
coomassie blue stain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors