

Recombinant Human SNAP23 protein ab79180

[1 References](#) [2 图像](#)

描述

产品名称	重组人SNAP23蛋白
纯度	> 90 % SDS-PAGE. ab79180 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MDNLSSEIIQ QRAHQITDES LESTRILGL AIESQDAGIK TITMLDEQKE QLNRIEEGLD QINKDMRETE KLTTELNKCC GLCVCPCNRT KNFESGKAYK TTWGDGGENS PCNVVSKQPG PVTNGQLQQP TTGAASGGYI KRITNDARED EMEENLTQVG SILGNLKDMA LNIGNEIDAQ NPQIKRITDK ADTNRDRIDI ANARAKKLID S

技术指标

Our **Abpromise guarantee** covers the use of **ab79180** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Western blot
形式	Liquid

制备和贮存


稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 8.00 Constituent: 0.242% Tris
--------	---

功能	Essential component of the high affinity receptor for the general membrane fusion machinery and an important regulator of transport vesicle docking and fusion.
组织特异性	Ubiquitous. Highest levels where found in placenta.
序列相似性	Belongs to the SNAP-25 family. Contains 2 t-SNARE coiled-coil homology domains.
细胞定位	Cell membrane. Cell membrane. Cell junction, synapse, synaptosome. Mainly localized to the plasma membrane.

图片



15% SDS-PAGE showing ab79180 (3µg).



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors