

Recombinant Human SQSTM1 / p62 protein ab95320

4 References 1 图像

|           |  |
|-----------|--|
| 描述        |  |
| 产品名称      | 重组人SQSTM1 / p62蛋白  |
| 纯度        | > 85 % SDS-PAGE.<br>ab95320 is purified using conventional chromatography techniques   |
| 表达系统      | Escherichia coli   |
| Accession | <u>Q13501</u>  |
| 蛋白长度      | Protein fragment   |
| 无动物成分     | No   |
| 性质        | Recombinant  |
| 种属        | Human  |
| 序列        | MAMSYVKDDI FRIYIKEKKE CRRDHRPPCA<br>QEAPRNMVHP NVICDGCNGP VVGTRYKCSV<br>CPDYDLCSVC EGKGLHRGHT KLAFPSPFGH<br>LSEGFSHSRW LRKVKHGHFG WPGWEMGPPG<br>NWSRPPRAG EARPGPТАES ASGPSEDPSV<br>NFLKNVGESV AAALSPLGIE VDIVVEHGGK<br>RSRLTPVSPE SSSTEEKSSS QPSSCCSDPS<br>KPGGNVEGAT QSLAEQMRKI ALESEGRPEE<br>QMESDNCSGG DDDWTHLSSK EVD PSTGELQ<br>SLQMPSEGP SSLDPSQEGP TGLKEAALYP<br>HLPPEADPRL IESLSQMLSM GFSDEGGWLT<br>RLLQTKNYDI GAALDTIQYS KHPPPLLEHH HHHH |
| 氨基酸       | 85 to 440  |
| 标签        | His tag C-Terminus   |

|  |                               |
|--|-------------------------------|
| 技术指标   |                               |
| Our <b>Abpromise guarantee</b> covers the use of <b>ab95320</b> in the following tested applications.                                |                               |
| The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user. |                               |
| 应用   | SDS-PAGE<br>Mass Spectrometry |
| 形式   | Liquid                        |

## 制备和贮存

---

### 稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

## 常规信息

---

### 功能

Adapter protein which binds ubiquitin and may regulate the activation of NFκB1 by TNF-α, nerve growth factor (NGF) and interleukin-1. May play a role in titin/TTN downstream signaling in muscle cells. May regulate signaling cascades through ubiquitination. Adapter that mediates the interaction between TRAF6 and CYLD (By similarity). May be involved in cell differentiation, apoptosis, immune response and regulation of K(+) channels.

### 组织特异性

Ubiquitously expressed.

### 疾病相关

Defects in SQSTM1 are a cause of Paget disease of bone (PDB) [MIM:602080]. PDB is a metabolic bone disease affecting the axial skeleton and characterized by focal areas of increased and disorganized bone turn-over due to activated osteoclasts. Manifestations of the disease include bone pain, deformity, pathological fractures, deafness, neurological complications and increased risk of osteosarcoma. PDB is a chronic disease affecting 2 to 3% of the population above the age of 40 years.

### 序列相似性

Contains 1 OPR domain.

Contains 1 UBA domain.

Contains 1 ZZ-type zinc finger.

### 结构域

The UBA domain binds specifically 'Lys-63'-linked polyubiquitin chains of polyubiquitinated substrates. Mediates the interaction with TRIM55.

The OPR domain mediates homooligomerization and interactions with PRKCZ, PRKCI, MAP2K5 and NBR1.

The ZZ-type zinc finger mediates the interaction with RIPK1.

### 翻译后修饰


Phosphorylated. May be phosphorylated by PRKCZ (By similarity). Phosphorylated in vitro by TTN.

### 细胞定位

Cytoplasm. Late endosome. Nucleus. Sarcomere (By similarity). In cardiac muscles localizes to the sarcomeric band (By similarity). Localizes to late endosomes. May also localize to the nucleus. Accumulates in neurofibrillary tangles and in Lewy bodies of neurons from individuals with Alzheimer and Parkinson disease respectively. Enriched in Rosenthal fibers of pilocytic astrocytoma. In liver cells, accumulates in Mallory bodies associated with alcoholic hepatitis, Wilson disease, indian childhood cirrhosis and in hyaline bodies associated with hepatocellular carcinoma.

## 图片

---



SDS-PAGE - Recombinant Human SQSTM1 / p62 protein (ab95320)

ab95320 at 3 µg on an SDS-PAGE gel (15%).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors