

Recombinant human TAK1 + TAB1 protein ab89692

5 图像

描述

产品名称	重组人TAK1 + TAB1蛋白
生物活性	The Specific activity of ab89692 was determined to be 10 nmol/min/mg.
纯度	> 80 % Densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human

技术指标

Our **Abpromise guarantee** covers the use of **ab89692** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Western blot Functional Studies
形式	Liquid
补充说明	ab64311 (Myelin Basic Protein protein) can be utilized as a substrate for assessing kinase activity


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.0038% EGTA, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	--

常规信息


功能	May be an important signaling intermediate between TGFB receptors and MAP3K7/TAK1. May play an important role in mammalian embryogenesis.
组织特异性	Ubiquitous.
序列相似性	Contains 1 PP2C-like domain.
翻译后修饰	Monoubiquitinated. Deubiquitinated by Y.enterocolitica YopP.

图片




The specific activity of TAK1 + TAB1 (ab89692) was determined to be 12 nmol/min/mg as per activity assay protocol

Functional Studies - Recombinant human TAK1 + TAB1 protein (ab89692)




SDS PAGE analysis of ab89692

SDS-PAGE - Recombinant human TAK1 + TAB1 protein (ab89692)




SDS PAGE analysis of ab89692

SDS-PAGE - Recombinant human TAK1 + TAB1
protein (ab89692)



SDS-PAGE showing ab89692 at approximately 74kDa.

SDS-PAGE - Recombinant human TAK1 + TAB1
protein (ab89692)



Kinase Assay demonstrating specific activity of ab89692.

Functional Studies - Recombinant human TAK1 +
TAB1 protein (ab89692)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors