

Recombinant Human TNFRSF19L protein ab191989

描述

产品名称	重组人TNFRSF19L蛋白	
纯度	> 95 % SDS-PAGE. determined by SEC-HPLC and reducing SDS-PAGE.	
内毒素水平	< 1.000 Eu/μg	
表达系统	HEK 293 cells	
Accession	<u>Q969Z4</u>	
蛋白长度	Protein fragment	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列		STTLWQCPPGEEPDLDPGQGTLCRPCPPGTFSAAWGSSPCQP HARCSLWR RLEAQVGMATRDTLGDCWPGWFGVPVPCQPCSWAPLG THGCDEWG RRARRGVEVAAGASSGGETRQPGNGTRAGGPEETA VDDIEGR MDEPKSCD KTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVV VDVSHEDP EVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDW LNGKEYKC KVS NKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQV SLTCLVKG FYPSDIAVEWESNGQPENNYKTTTPVLDSGDGFFLYSKLTVD KSRWQQGN VFSCVMHEALHNHYTQKSLSLSPGK
预测分子量	41 kDa including tags	
氨基酸	26 to 160	
标签	Fc tag C-Terminus	
额外的序列信息	Extracellular domain, fused with a FC tag at the C-terminus.	

技术指标

Our **Abpromise guarantee** covers the use of **ab191989** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	HPLC SDS-PAGE
形式	Liquid

制备和贮存

稳定性和存储	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. pH: 7.2 Constituents: 1.19% HEPES, 0.87% Sodium chloride, 10% Glycerol (glycerin, glycerine), 0.03% DTT
--------	---

常规信息

功能	Mediates activation of NF-kappa-B. May play a role in T-cell activation.
组织特异性	Highest levels are in spleen, lymph node, thymus, peripheral blood leukocytes, bone marrow and fetal liver. Very low levels in skeletal muscle, testis and colon. Not detected in brain, kidney and pancreas.
序列相似性	Belongs to the RELT family. Contains 1 TNFR-Cys repeat.
翻译后修饰	Phosphorylated in vitro by OXSR1.
细胞定位	Cell membrane. Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors