

Recombinant Human TOMM20 protein ab156966

1 图像

描述	
产品名称	重组人TOMM20蛋白
纯度	> 90 % SDS-PAGE. ab156966 was purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q15388</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMGSDRKRRSDPNFKNRLRERRK KQKLAKER AGLSKLPDLKDAAEQKFFLEEIQLGEELLAQGEYEKGV DHL TNAIAVCG QPQQLLQVLQQTLP PPVFQMLLTKLPTISQRIVSAQSLAEDD VE
预测分子量	16 kDa including tags
氨基酸	25 to 145
标签	His tag N-Terminus
技术指标	
Our Abpromise guarantee covers the use of ab156966 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Mass Spectrometry SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.04% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

功能	Central component of the receptor complex responsible for the recognition and translocation of cytosolically synthesized mitochondrial preproteins. Together with TOM22 functions as the transit peptide receptor at the surface of the mitochondrion outer membrane and facilitates the movement of preproteins into the TOM40 translocation pore.
序列相似性	Belongs to the Tom20 family.
细胞定位	Mitochondrion outer membrane.

图片



15% SDS-PAGE analysis of ab156966 (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors