

Recombinant Human Ubiquilin/UBQLN1 protein ab153310

1 图像

描述	
产品名称	重组人Ubiquilin/UBQLN1蛋白
表达系统	Wheat germ
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAESGESGGPPGSQDSAAGAEGAGTPAAAASAEPKIMKVTVK TPKEKEEF AVPENSSVQQFKEEISKRFKSHTDQLVLIFAGKILKDQDTLS QHGIHDGL TVHLVIKTQNRPDHSAQQTNTAGSNVTTSSTPNSNSTSGSA TSNPFGLG GLGGLAGLSSLGLNTTNFSELQSQMQRQLLSNPEMMVQIMEN PFVQSMLS NHDLMRQLIMANPQMQLIQRNPEISHMLNPNPDIMRQTLELA RNPAMMQE MMRNQDRALSNLESIPGGYNALRRMYTDIQEPMLSAAQEQFG GNPFASLV SNTSSGEGSQPSRTENRDPLPNPWAPQTSQSSSASSGTASTV GGTTGSTA SGTSGQSTTAPNLVPGVGASMFNTPGMQSLLQQITENPQLMQ NMLSAPYM RSMMQSLSQNPDLAAQMMLNNPLFAGNPQLQEQMRQQLPTFL QQMQNPDT LSAMSNPRAMQALLQIQQGLQTLATEAPGLIPGFTPGLGALG STGGSSGT NGSNATPSENTSPTAGTTEPGHQQFIQQMLQALAGVNPQLQN PEVRFQQQ LEQLSAMGFLNREANLQALIATGGDINAAIERLLGSQPS
氨基酸	1 to 589
标签	GST tag N-Terminus
技术指标	

Our **Abpromise guarantee** covers the use of **ab153310** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot ELISA
形式	Liquid
补充说明	This product was previously labelled as Ubiquilin.


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
--------	--

常规信息

功能	Links CD47 to the cytoskeleton. Promotes the surface expression of GABA-A receptors (By similarity). Promotes the accumulation of uncleaved PSEN1 and PSEN2 by stimulating their biosynthesis. Has no effect on PSEN1 and PSEN2 degradation.
组织特异性	Ubiquitous. Highly expressed throughout the brain; detected in neurons and in neuropathological lesions, such as neurofibrillary tangles and Lewy bodies. Highly expressed in heart, placenta, pancreas, lung, liver, skeletal muscle and kidney.
序列相似性	Contains 1 UBA domain. Contains 1 ubiquitin-like domain.
结构域	Binds PSEN1 and PSEN2 via its UBA domain.
翻译后修饰	Phosphorylated.
细胞定位	Cytoplasm. Nucleus.

图片



ab153310 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors