

Recombinant Human VPAC2 protein ab159805

1 图像

描述	
产品名称	重组人VPAC2蛋白
表达系统	Wheat germ
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MRTLLPPALLTCWLLAPVNSIHPECRFHLEIQEEETKCAELL RSQTEKHK ACSGVWDNITCWRPANVGETVTVPCPKVFSNFYSKAGNISK CTSDGWSE TFPDFVDACGYSDPEDESKITFYILVKAIYTLGYSVSLMSLA TGSIIILCL FRKLHCTRNYIHLNLFSLFILRAISVLVKDDVLYSSSGTLHC PDQPSSWV GCKLSLVFLQYCIMANFFWLLVEGLYLHTLLVAMLP PRRCFLAYLLIGWG LPTVCIGAWTAARLYLEDTGCDTNDHSVPWWVIRIPILISI IVNFVLF SIIRILLQKLTPDVGGNDQSQYKRLAKSTLLLIPLFGVHYM VFAVFPIS ISSKYQILFELCLGSFQGLVVAVLYCFLNSEVQCE LKRKWSRCPTPSAS RDYRVCGSSFSRNGSEGALQFHRGSRAQSFLQTETSVI
氨基酸	1 to 438
标签	GST tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab159805 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	ELISA Western blot

形式 Liquid

补充说明

制备和贮存

稳定性和存储 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

常规信息


功能 This is a receptor for VIP as well as PACAP-38 and -27, the activity of this receptor is mediated by G proteins which activate adenylyl cyclase. Can be coupled to phospholipase C.

组织特异性 Expressed in CD4+ T-cells, but not in CD8+ T-cells. Expressed in the T-cell lines Jurkat, PEER, MOLT-4, HSB, YT and Tsup-1, but not in the T-cell lines HARRIS and HUT 78.

序列相似性 Belongs to the G-protein coupled receptor 2 family.

细胞定位 Cell membrane.

图片



SDS-PAGE - Recombinant Human VPAC2 protein
(ab159805)

ab159805 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors