

Recombinant Mouse CNDP1/CN1 protein (His tag) ab222350

1 图像

描述

产品名称	重组小鼠CNDP1/CN1蛋白(His tag)
纯度	> 95 % SDS-PAGE. Affinity purified
内毒素水平	< 1.000 Eu/μg
表达系统	Baculovirus infected insect cells
Accession	Q8BUG2
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse
序列	<p>MFSSAHSGLLEKLFHYIDLHQDEFVQTLKEWVAIESDSVQPV PRLRQKLF QMMALAADKLRNLGAGVESIDLGSQQMPDGQSLPIPPILLAE LGSDPEKP TVCFYGHLDVQPAQKDDGWLTPYTLTEVDGKLYGRGATDNK GPVLAWIN AVSTFRALQQDLPVNIKLI LEGMEEAGSIALEELVMREKDH FSSVDYIV ISDNLWLSQRKPALTYGTRGNCYFTVEVKCRDQDFHSGTFGG ILNEPMAD LVALLGSLVDSSGHILIPGIYDQMAPITEGEKTM YKNIDMDL EEYQNIQ VEKFLFDTKEELLMHLWRYPSLSIHGIEGAFDEPGTKTVIPG RVLGKFSI RLVPTMSPSVVEKQVTQHLEAVFSKRNSFNKMAVSMVLGLHP WTANVNDT QYLAAQRTIKTVFGVNPDMIRDGSTIPIAKIFQAITQKSVMM LPLGAVDD GEHSQNEKINRWNYIQGSKLFAAFFLELSKQHSQHGHPSSVY LEHHHHHH</p>
预测分子量	56 kDa including tags
氨基酸	1 to 492
标签	His tag C-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab222350** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid
补充说明	This product was previously labelled as CNDP1

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.40 Constituents: 10% Glycerol (glycerin, glycerine), PBS
--------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

常规信息

组织特异性	Found in serum and adult nervous central system. Absent in serum from patients with homocarnosinosis.
序列相似性	Belongs to the peptidase M20A family.
细胞定位	Secreted.

图片



15% SDS-PAGE analysis of 3 µg ab222350.

MWt 50-70 kDa (SDS-PAGE under reducing conditions).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors