

Anti-TMPRSS2 antibody [EPR3861] Blocking peptide ab275042

1 图像

描述

产品名称	Anti-TMPRSS2抗体[EPR3861] Blocking多肽
纯度	> 90 % HPLC.
Accession	<u>O15393</u>
无动物成分	No
性质	Synthetic
种属	Human
描述	人TMPRSS2多肽

技术指标

Our **Abpromise guarantee** covers the use of **ab275042** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking
形式	Lyophilized

制备和贮存





稳定性和存储	Shipped at 4°C. Store at -20°C.
复溶	Reconstitute in either water or buffer. If the peptide doesn't dissolve try an organic solvent like DMSO, then dilute using water or buffer. Gentle warming and sonication can effectively aid peptide solubilization.

常规信息

组织特异性	Expressed strongly in small intestine. Also expressed in prostate, colon, stomach and salivary gland.
序列相似性	Belongs to the peptidase S1 family. Contains 1 LDL-receptor class A domain. Contains 1 peptidase S1 domain. Contains 1 SRCR domain.

图片

Choose the best for your experiment from thousands of proteins and peptides

-  Reinforce the validity of your results with **protein controls**
-  Validate specific, reliable reagents with **blocking peptides**
-  Keep working on targets without specific antibodies with **fusion-tagged proteins**
-  Get consistent, reproducible results with premium **bioactive proteins**

Anti-TMPRSS2 antibody [EPR3861] Blocking peptide (ab275042)

To learn more about our protein and peptide range click [here](#).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors