

## Product datasheet

### Guinea Pig Complement Preserved Serum ab155111

#### 概述

---

产品名称	豚鼠Complement Preserved Serum
宿主	Guinea pig
常规说明	Guinea Pig Complement serum is collected from spontaneously clotted whole blood and centrifuged to separate the serum from the clot. Guinea Pig Complement is ideal for hemolytic assay where hemolysis of IgG antibodies are required. Complement is preservative free. Raised in: Guinea pig Purity: Whole antiserum

#### 性能

---

形式	Liquid
存放说明	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors