



Product datasheet

Normal Hamster Serum (sterile) ab7483

概述

| | |
|-------|---|
| 产品名称 | Normal仓鼠Serum (sterile) |
| 宿主 | Syrian hamster |
| 经测试应用 | 适用于: Blocking, IHC-P, IHC-Fr |
| 常规说明 | Total Protein Concentration : normal Strain: Golden Syrian Sex: Mixed sex Raised in: Syrian Hamster Purity: Whole antiserum |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7483于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----------|------|--|
| Blocking | | Use at an assay dependent concentration. |
| IHC-P | | Use at an assay dependent concentration. |
| IHC-Fr | | Use at an assay dependent concentration. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors