

Donkey Anti-Goat IgG H&L (Cy3 ®) preadsorbed ab6949

★★★★★ [2 Abreviews](#) [41 References](#) [1 图像](#)

概述

产品名称	驴抗山羊IgG H&L (Cy3 ®)预吸附二抗
宿主	Donkey
靶标种属	Goat
经测试应用	适用于: ICC/IF, IHC-P, IHC-Fr, ICC, Flow Cyt, ELISA, WB
最小交叉反应	Cow, Guinea pig, Hamster, Horse, Mouse, Rabbit, Rat more details
免疫原	Full length native Goat IgG (purified).
偶联物	Cy3 ®. Ex: 552nm, Em: 565nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA
纯度	Immunogen affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities and extensive dialysis against the buffer stated above.
共轭说明	Cy3.29 (Cyanine 3.29-OSu) (Molecular Weight 949 daltons) Absorption Wavelength: 552 nm Emission Wavelength: 565 nm Fluorochrome/Protein Ratio: 8.1 moles Cy3 per mole of Donkey IgG
克隆	多克隆
同种型	IgG
常规说明	Cy™ and CyDye™ are registered trademarks of Cytiva. This secondary antibody is specifically designed for the detection of multiple primary antibodies (polyclonal or monoclonal) of different host species in experiments where cells are simultaneously labeled without unwanted cross reaction.


应用

The Abpromise guarantee Abpromise™ 承诺保证使用ab6949于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/1000 - 1/5000.
IHC-P		Use at an assay dependent concentration.
IHC-Fr	★★★★★ (1)	Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration.
ELISA		1/10000 - 1/50000.
WB		1/1000. Predicted molecular weight: 36 kDa.

图片



Donkey Anti-Goat IgG H&L (Cy3 ®) preadsorbed (ab6949) at 1000 cells + Goat IgG, 50 ng

Predicted band size: 36 kDa

Both the heavy and light chains of Goat IgG are observed.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors