

Donkey Anti-Goat IgG H&L (Texas Red®) preadsorbed ab7123

2 References

概述

产品名称	驴抗山羊IgG H&L (Texas Red®)预吸附二抗
宿主	Donkey
靶标种属	Goat
经测试应用	适用于: ICC/IF, Immunomicroscopy, Flow Cyt, ELISA, IHC-P, IHC-Fr
最小交叉反应	Chicken, Guinea pig, Hamster, Horse, Mouse, Rabbit, Rat more details
免疫原	Full length native Goat IgG (purified).
偶联物	Texas Red®. Ex: 596nm, Em: 620nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Texas Red™ Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 2.1 moles Texas Red® per mole of Rabbit IgG
克隆	多克隆

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab7123于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		
Immunomicroscopy		
Flow Cyt		
ELISA		
IHC-P		
IHC-Fr		

应用说明

Flow Cyt: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IM: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors