



Product datasheet

Donkey Anti-Mouse IgG H&L preadsorbed ab7056

2 References

概述		
产品名称	驴抗小鼠IgG H&L预吸附二抗	
宿主	Donkey	
靶标种属	Mouse	
经测试应用	适用于: WB, IP, ICC/IF, Immunodiffusion, Conjugation, IHC-P, IHC-Fr, ELISA	
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep	
免疫原	Full length native Mouse IgG (purified).	more details

性能		
形式	Liquid	
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.	
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride	
纯度	Affinity purified	
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose	
克隆	多克隆	

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab7056于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/2000 - 1/10000.
IP		Use at an assay dependent concentration.

应用	Ab评论	说明
ICC/IF		1/1000 - 1/5000.
Immunodiffusion		Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ELISA		1/30000 - 1/120000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors