

Donkey Anti-Rabbit IgG H&L (TRITC) preadsorbed ab7080

[1 References](#) [1 图像](#)

概述

产品名称	驴抗兔IgG H&L (TRITC)预吸附二抗
宿主	Donkey
靶标种属	Rabbit
经测试应用	适用于: ICC/IF, Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr, ELISA
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rat, Sheep more details
免疫原	Full length native Rabbit IgG (purified).
偶联物	TRITC. Ex: 547nm, Em: 572nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	Tetramethylrhodamine isothiocyanate (TRITC) (Molecular Weight 444 daltons) Absorption Wavelength: 550 nm Emission Wavelength: 570 nm Fluorochrome/Protein Ratio: 2.3 moles TRITC per mole of Goat IgG
共轭说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7080于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		
Immunomicroscopy		
Flow Cyt		
IHC-P		
IHC-Fr		
ELISA		

应用说明

Flow Cyt: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.


IHC-Fr: Use at an assay dependant concentration.

IM: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

图片



ab7080 was used at dilution 1/200 with the primary antibody **ab2900** in ICC/IF. See the review on **ab2900**.

Immunocytochemistry/ Immunofluorescence -
Donkey Anti-Rabbit IgG H&L (TRITC) preadsorbed
(ab7080)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors