



Product datasheet

Goat Anti-Chicken IgY H&L (25nm Gold) ab41494

概述

产品名称	山羊抗鸡IgY H&L (25nm Gold)
宿主	Goat
靶标种属	Chicken
经测试应用	适用于: Electron Microscopy
免疫原	Chicken IgG whole molecule.
偶联物	Gold 25nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.097% Sodium azide Constituents: PBS, 1% BSA
纯度	Immunogen affinity purified
纯化说明	This product has been immuno-affinity purified and immuno cross-absorbed to reduce non-specific reactions. The activity of each lot is determined using a dot-spot test system [Moeremans et al., J. Immunol., Methods 74, (1984), 353.]
共轭说明	This product is an immunogold reagent built around colloidal gold particles of 6nm. The particle population is monodisperse and thus shows minimal size variation and overlap. Typically, the coefficient of variance for the 6nm particle size conjugate is less than 12%.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee Abpromise™承诺保证使用ab41494于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Electron Microscopy		

应用说明

Electron Microscopy: Use at an assay dependent dilution. Electron Microscopy: Use at an assay dependent dilution. The method of production of this product allows the small particle size of 6nm to be used instead of 5nm with the same sensitivity but increased visibility in the standard transmission electron microscope (See also van Bergen en Henegouwen and Leunissen, Histochemistry, 85,(1986), 81.)

Light Microscopy: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors