



Product datasheet

Goat Anti-Chicken IgY H&L preadsorbed ab7113

1 References

概述		
产品名称	山羊抗鸡IgY H&L预吸附二抗	
宿主	Goat	
靶标种属	Chicken	
经测试应用	适用于: WB, ICC/IF, IHC-Fr, IHC-P, IP, Immunodiffusion, Conjugation, ELISA	
最小交叉反应	Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rabbit, Rat, Sheep	
免疫原	Full length native Chicken IgY (purified).	
		<a href="#">more details</a>

性能		
形式	Liquid	
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.	
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride	
纯度	Affinity purified	
纯化说明	Prepared from monospecific antiserum by immunoaffinity chromatography using Chicken IgG coupled to agarose beads.	
克隆	多克隆	

应用

The Abpromise guarantee      **Abpromise™** 承诺保证使用ab7113于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/2000 - 1/10000.
ICC/IF		1/1000 - 1/5000.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.
ELISA		1/600000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors