

Goat Anti-Dog IgG H&L ab112849

概述	
产品名称	山羊抗狗IgG H&L
宿主	Goat
靶标种属	Dog
特异性	This antibody reacts specifically with dog IgG and with light chains common to other dog immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species
经测试应用	适用于: WB, ICC, IHC, ELISA
免疫原	Full length native protein (purified) corresponding to Dog IgG. Purified dog IgG

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab112849 was isolated by affinity chromatography using antigen coupled to agarose beads.
克隆	多克隆
同种型	IgG

应用		
The Abpromise guarantee <u>Abpromise™</u> 承诺保证使用ab112849于以下的经测试应用		
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。		
应用	Ab评论	说明
WB		1/1000 - 1/3000.

应用	Ab评论	说明
ICC		1/200 - 1/2000.
IHC		1/200 - 1/2000.
ELISA		1/1000 - 1/30000. Coating plates: use 1:100 - 1:500.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors