

Goat F(ab) Anti-Mouse IgG H&L ab6668

36 References 2 图像

概述

产品名称	山羊F(ab)抗小鼠IgG H&L
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: WB, ICC/IF, IHC, ELISA
免疫原	Mouse IgG whole molecule

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, papain digestion and chromatographic separation. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum. No reaction was observed against anti-Papain or anti-Goat IgG F(c).”
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee Abpromise™承诺保证使用ab6668于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/2000 - 1/10000.
ICC/IF		1/1000 - 1/5000.

应用	Ab评论	说明
IHC		1/1000 - 1/5000.
ELISA		1/50000.

图片



SDS-PAGE results of Goat Fab Anti-Mouse IgG (H&L) ab6668.


Lane 1: Non-reduced Goat Fab Anti-Mouse.

Lane 2: Opal Prestained Ladder.

Lane 3: Reduced Goat Fab Anti-Mouse.

Load: 1.0 µg.

Coomassie stained and imaged.



Dot Blot results of Goat Fab Anti-Mouse IgG Antibody ab6668.

Dots are Mouse IgG at (1) 100ng, (2) 33.3ng, (3) 11.1ng, (4) 3.70ng, (5) 1.23ng.

Blocked for 30 min at RT.

Primary Antibody: Goat Fab Anti-Mouse IgG at 1 µg/mL for 1 hr at RT.

Secondary Antibody: Donkey Anti-Goat IgG HRP at 1:40,000 for 30 min at RT.

Imaged with BioRad ChemiDoc, Chemi filter.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors