

Goat F(ab) Anti-Rabbit IgG H&L (Biotin) ab7055

1 References **1 图像**

概述

| | |
|-------|---|
| 产品名称 | 山羊F(ab)抗兔IgG H&L (Biotin) |
| 宿主 | Goat |
| 靶标种属 | Rabbit |
| 特异性 | Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Biotin and anti-Goat Serum. No reaction was observed against anti-Papain or anti-Goat IgG F(c). |
| 经测试应用 | 适用于: ICC/IF, Dot blot, WB, ELISA, IHC, Immunomicroscopy |
| 免疫原 | Rabbit IgG whole molecule |
| 偶联物 | Biotin |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA |
| 纯化说明 | BSA is IgG and protease free This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, papain digestion and chromatographic separation. |
| 共轭说明 | Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Goat IgG Fab molecule |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab7055于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|------------------|------|----|
| ICC/IF | | |
| Dot blot | | |
| WB | | |
| ELISA | | |
| IHC | | |
| Immunomicroscopy | | |

应用说明

Dot: Use at an assay dependent dilution.

ELISA: 1/50000 - 1/250000.

IHC: Use at an assay dependent dilution.

IM: Use at an assay dependent dilution.

WB: Use at an assay dependent dilution.


This product has been assayed against 1.0 µg of Rabbit IgG in a standard capture ELISA using Peroxidase Conjugated Streptavidin and ABTS (2,2'-azino-bis-[3-ethylbenzothiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature. A working dilution of 1/50000 to 1/250000 of the reconstitution concentration is suggested for this product.

Suitable for other antibody based assays using streptavidin or avidin conjugates requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat F(ab) Anti-Rabbit IgG H&L (Biotin) (ab7055)

ab7055 was used at dilution 1/200 with the primary antibody **ab32090** in IHC-P. See the review on **ab32090**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors