

## Product datasheet

# Goat F(ab')<sub>2</sub> Anti-Human IgG - Fc (DyLight® 550), pre-adsorbed ab98591

### 1 图像

#### 概述

产品名称	山羊F(ab') <sub>2</sub> 抗人IgG - Fc (DyLight® 550),预吸附二抗
宿主	Goat
靶标种属	Human
经测试应用	适用于: Flow Cyt, IHC-P, ICC/IF
最小交叉反应	Mouse, Rat
偶联物	DyLight® 550. Ex: 562nm, Em: 576nm

[more details](#)

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 550.
克隆	多克隆
同种型	IgG
常规说明	DyLight® 550 replaces DyLight® 549.

#### 应用


##### The Abpromise guarantee

**Abpromise™** 承诺保证使用ab98591于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		1/50 - 1/200.
IHC-P		1/50 - 1/500.
ICC/IF		1/50 - 1/500.

## 图片



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors