

Goat F(ab')₂ Anti-Human IgG - Fc (DyLight® 594), pre-adsorbed ab98592

1 图像


概述	
产品名称	山羊F(ab') ₂ 抗人IgG - Fc (DyLight® 594),预吸附二抗
宿主	Goat
靶标种属	Human
经测试应用	适用于: Flow Cyt, ICC/IF, IHC-P
最小交叉反应	Mouse, Rat
偶联物	DyLight® 594. Ex: 593nm, Em: 618nm
more details	

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 594.
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	<u>Abpromise™</u> 承诺保证使用ab98592于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
Flow Cyt		1/50 - 1/200.
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.

图片



Goat F(ab')₂ Anti-Human IgG - Fc (DyLight® 594),
pre-adsorbed (ab98592)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors