

# Goat F(ab')<sub>2</sub> Anti-Human IgG H&L (PE) preadsorbed ab7006

## 6 References

### 概述

产品名称	山羊F(ab') <sub>2</sub> 抗人IgG H&L (PE)预吸附二抗	
宿主	Goat	
靶标种属	Human	
经测试应用	适用于: ICC/IF, IHC-P, IHC-Fr, Immunomicroscopy, Flow Cyt	
最小交叉反应	Cow, Horse, Mouse, Rat	<a href="#">more details</a>
免疫原	Full length native Human IgG (purified).	
偶联物	PE. Ex: 488nm, Em: 575nm	

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.05% Sodium azide Constituents: 0.88% Sodium chloride, 1% BSA, 0.42% Potassium phosphate
纯度	Immunogen affinity purified
纯化说明	Prepared from monospecific antiserum by immunoaffinity chromatography using human IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation
共轭说明	R-Phycoerythrin (R-PE) Molecular Weight 240,000 daltons) Absorption Wavelength: 488 nm Emission Wavelength: 575 nm OD <sub>566nm</sub> /OD <sub>280nm</sub> : 2.8 Phycobiliprotein Concentration: 0.5 mg / ml (determined by absorbance = 82.0 at 565 nm for a 1% solution for only those R-PE molecules to which at least one molecule of active antibody is bound).
克隆	多克隆
同种型	IgG
常规说明	Can be used in immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.

## The Abpromise guarantee **Abpromise™**承诺保证使用ab7006于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/100 - 1/250.
IHC-P		1/50 - 1/200.
IHC-Fr		1/50 - 1/200.
Immunomicroscopy		1/50 - 1/200.
Flow Cyt		1/100 - 1/250. Use $1 \times 10^6$ cells in flow cytometry is approximately 1.0 ug of antibody conjugate.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors