

# Goat F(ab')<sub>2</sub> Anti-Mouse IgG - F(ab')<sub>2</sub> (HRP), pre-adsorbed ab5887

## 7 References

### 概述

产品名称	山羊F(ab') <sub>2</sub> 抗小鼠IgG - F(ab') <sub>2</sub> (HRP),预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat Serum, Mouse IgG, Mouse IgG F(ab') <sub>2</sub> and Mouse Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c), Mouse IgG F(c) or Bovine, Horse, Human, Rabbit, Rat and Sheep Serum Proteins.
经测试应用	适用于: ICC/IF, Dot blot, WB, ELISA, IHC, Immunomicroscopy
最小交叉反应	Cow, Horse, Human, Rabbit, Rat, Sheep
免疫原	Mouse IgG F(ab) <sub>2</sub> fragment
偶联物	HRP

[more details](#)

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯化说明	BSA is IgG and protease free  This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
共轭说明	Horseradish Peroxidase (HRP) (Molecular Weight 40,000 daltons)
克隆	多克隆
同种型	IgG

## 应用

**The Abpromise guarantee** Abpromise™承诺保证使用ab5887于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/500 - 1/2500.
Dot blot		Use at an assay dependent dilution.
WB		1/2000 - 1/20000.
ELISA		1/100000. This product has been assayed against 1.0 ug of Mouse IgG in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature.
IHC		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors