

Goat F(ab')₂ Anti-Mouse IgG+IgM+IgA - H&L (AP), pre-adsorbed ab6007

概述

产品名称	山羊F(ab') ₂ 抗小鼠IgG+IgM+IgA - H&L (AP),预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase and anti-Goat Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c) or Human Serum Proteins. This reagent is suitable for the detection of all mouse immunoglobulin subclasses and chain combinations.
经测试应用	适用于: ICC/IF, Dot blot, WB, ELISA, IHC
最小交叉反应	Human
免疫原	Mouse IgG, IgA and IgM whole molecules
偶联物	Alkaline Phosphatase

[more details](#)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 1% BSA
纯化说明	BSA is IgG and protease free This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons).
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee Abpromise™承诺保证使用ab6007于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/100 - 1/500.
Dot blot		Use at an assay dependent dilution.
WB		1/200 - 1/1000.
ELISA		1/6300.
IHC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors