

Goat F(ab')₂ Anti-Mouse IgG+IgM+IgA H&L (APC) preadsorbed ab130775

概述

产品名称	山羊F(ab') ₂ 抗小鼠IgG+IgM+IgA H&L (APC)预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	ab130775 reacts with the heavy and light chains of mouse IgM, IgG1, IgG2a, IgG2b IgG3, and IgA. Goat F(ab') ₂ Anti-Mouse IgG+IgM+IgA H&L (Allophycocyanin) preadsorbed shows minimal cross reactivity with Human immunoglobulins.
经测试应用	适用于: ICC/IF, ELISA, WB, Dot blot, IHC-P, IHC-Fr
最小交叉反应	Human more details
免疫原	Full length native protein (purified) corresponding to Mouse IgG + IgM + IgA. Mouse Ig
偶联物	APC. Ex: 645nm, Em: 660nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.3 Preservative: 0.1% Sodium azide Constituent: PBS
纯化说明	Sucrose is added as a stabilizing agent. ab130775 was purified by gel filtration chromatography of pepsin digested antibody.
克隆	多克隆
同种型	IgG
常规说明	ab130775 is adsorbed against pooled Human sera and purified Human paraproteins.

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab130775于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use 0.1µg for 10 ⁶ cells.
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors