

Goat F(ab')₂ Anti-Mouse IgM mu chain (HRP) ab5930

12 References

概述

产品名称	山羊F(ab') ₂ 抗小鼠IgM mu chain (HRP)
宿主	Goat
靶标种属	Mouse
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat Serum, Mouse IgM and Mouse Serum. No reaction was observed against anti-Pepsin or anti-Goat IgG F(c). Specificity was confirmed by ELISA at less than 1% cross reactivity against other mouse heavy or light chain isotypes.
经测试应用	适用于: ICC/IF, IHC-Fr, Dot blot, WB, ELISA, Electron Microscopy, IHC
免疫原	Mouse IgM mu chain
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.5 Preservative: 0.01% Gentamicin sulphate Constituents: 0.87% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgM coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
共轭说明	Horseradish Peroxidase (HRP) (Molecular Weight 40,000 daltons)
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab5930于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/500 - 1/5000.
IHC-Fr		Use at an assay dependent dilution.
Dot blot		1/5000 - 1/100000.
WB		1/2000 - 1/10000.
ELISA		1/90000.
Electron Microscopy		Use at an assay dependent dilution.
IHC		1/500 - 1/5000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors