

### Goat F(ab')<sub>2</sub> Anti-Rat IgG - (Fab')<sub>2</sub> (AP), pre-adsorbed ab6518

#### 概述

---

产品名称	山羊F(ab') <sub>2</sub> 抗大鼠IgG - (Fab') <sub>2</sub> (AP),预吸附二抗
宿主	Goat
靶标种属	Rat
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase, anti-Goat Serum, Rat IgG F(ab) <sub>2</sub> and Rat Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c), Rat Fc or Bovine, Horse, and Human Serum Proteins
经测试应用	适用于: ELISA, IHC, ICC/IF, Dot blot, WB
最小交叉反应	Cow, Horse, Human <a href="#">more details</a>
免疫原	Rat IgG F(ab) <sub>2</sub> fragment
偶联物	Alkaline Phosphatase

#### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium azide Constituents: 0.79% Tris HCl, 0.88% Sodium chloride, 1.5% BSA, 0.01% Magnesium chloride, 0.0014% Zinc chloride, 50% Glycerol (glycerin, glycerine)
纯化说明	BSA is IgG and protease free This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆
同种型	IgG

#### 应用

---

## The Abpromise guarantee

**Abpromise™**承诺保证使用ab6518于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/5000 - 1/50000.
IHC		Use at an assay dependent dilution.
ICC/IF		1/200 - 1/1000.
Dot blot		Use at an assay dependent dilution.
WB		1/5000 - 1/50000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors