

# Goat F(ab')<sub>2</sub> Anti-Rat IgG - (Fab')<sub>2</sub> (FITC), pre-adsorbed ab6266

## 2 References

### 概述

产品名称	山羊F(ab') <sub>2</sub> 抗大鼠IgG - (Fab') <sub>2</sub> (FITC),预吸附二抗
宿主	Goat
靶标种属	Rat
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein, anti-Goat Serum, Rat IgG, Rat IgG F(ab') <sub>2</sub> and Rat Serum. No reaction was observed against Rat IgG F(c), anti-Pepsin, anti-Goat IgG F(ab') <sub>2</sub> or Bovine, Horse and Human Serum Proteins.
经测试应用	适用于: ICC/IF, Flow Cyt, Immunomicroscopy, ELISA
最小交叉反应	Cow, Horse, Human <a href="#">more details</a>
免疫原	Rat IgG F(ab') <sub>2</sub> fragment
偶联物	FITC. Ex: 493nm, Em: 528nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.05% Sodium azide Constituents: 1.46% Sodium chloride, 1.5% BSA, 0.164% Sodium phosphate
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
共轭说明	Fluorescein isothiocyanate (FITC) (Molecular Weight 390 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 3.1 moles FITC per mole of Goat IgG F(ab') <sub>2</sub>
克隆	多克隆
同种型	IgG

### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab6266于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		
Flow Cyt		
Immunomicroscopy		
ELISA		

#### 应用说明

Flow Cyt: Use at an assay dependent dilution.

IM: Use at an assay dependent dilution.

Suitable for other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors